

SMARTINTERNZ

GUIDED PROJECT

ECG- Image Based Heartbeat Classification For
Arrhythmia Detection

Using IBM Watson Studio

- UPPANAPALLI LAKSHMI SOWJANYA

Project Description

According to the World Health Organization (WHO), cardiovascular diseases (CVDs) are the number one cause of death today. Over 17.7 million people died from CVDs in the year 2017 all over the world which is about 31% of all deaths, and over 75% of these deaths occur in low and middle-income countries. Arrhythmia is a representative type of CVD that refers to any irregular change from the normal heart rhythms. There are several types of arrhythmia including atrial fibrillation, premature contraction, ventricular fibrillation, and tachycardia. Although a single arrhythmia heartbeat may not have a serious impact on life, continuous arrhythmia beats can result in fatal circumstances. In this project, we build an effective electrocardiogram (**ECG**) arrhythmia classification method using a convolutional neural network (**CNN**), in which we classify ECG into seven categories, one being normal and the other six being different types of arrhythmia using deep two-dimensional CNN with **grayscale ECG images**. We are creating a web application where the user selects the image which is to be classified. The image is fed into the model that is trained and the cited class will be displayed on the webpage.

Project Objectives


By the end of this project I

- knew fundamental concepts and techniques of the Artificial Neural Network and Convolution Neural Networks
- Gained a broad understanding of image data.
- Worked with Sequential type of modeling
- Worked with Keras capabilities
- Worked with image processing techniques
- knew how to build a web application using the Flask framework.

Literature Review

Problem

Cardiovascular diseases (CVDs) are the number one cause of death today. Over 17.7 million people died from CVDs in the year 2017 all over the world which is about 31% of all deaths, and over 75% of these deaths occur in low and middle-income countries. Arrhythmia is a representative type of CVD that refers to any irregular change from the normal heart rhythms. There are several types of arrhythmia including atrial fibrillation, premature contraction, ventricular fibrillation, and tachycardia.




ARRHYTHMIA CAUSES

- ◆ Blocked arteries in the heart
- ◆ Scarring from a previous heart attack
- ◆ Changes to the heart's structure
- ◆ Diabetes
- ◆ High blood pressure
- ◆ Overactive thyroid gland
- ◆ Sleep apnea
- ◆ Underactive thyroid gland
- ◆ Certain medications
- ◆ Stress or anxiety
- ◆ Smoking

Solution

An "ambulatory electrocardiogram" or an ECG (about the size of a postcard or digital camera) that the patient will be using for 1 to 2 days, or up to 2 weeks. The test measures the movement of electrical signals or waves through the heart. These signals tell the heart to contract (squeeze) and pump blood. The patient will have electrodes taped to your skin. It's painless, although some people have mild skin irritation from the tape used to attach the electrodes to the chest. They can do everything but shower or bathe while wearing the electrodes. After the test period, patient will go back to see your doctor. They will be downloading the information.



Role of Deep Learning

Deep learning is a branch of machine learning which is completely based on artificial neural networks, as neural network is going to mimic the human brain so deep learning is also a kind of mimic of human brain. In deep learning, we don't need to explicitly program everything. The concept of deep learning is not new. It has been around for a couple of years now. It's on hype nowadays because earlier we did not have that much processing power and a lot of data. As in the last 20 years, the processing power increases exponentially, deep learning and machine learning came in the picture.


Deep learning is a particular kind of machine learning that achieves great power and flexibility by learning to represent the world as a nested hierarchy of concepts, with each concept defined in relation to simpler concepts, and more abstract representations computed in terms of less abstract ones.

CNN

In the field of deep learning, convolutional neural network (CNN) is among the class of deep neural networks, which was being mostly deployed in the field of analyzing/image recognition. Convolutional Neural uses a very special kind of method which is being known as Convolution.

The Convolutional neural networks(CNN) consists of various layers of artificial neurons. Artificial neurons, similar to that neuron cells that are being used by the human brain for passing various sensory input signals and other responses, are mathematical functions

that are being used for calculating the sum of various inputs and giving output in the form of an activation value.



The behaviour of each CNN neuron is being defined by the value of its weights. When being fed with the values (of the pixel), the artificial neurons of a CNN recognizes various visual features and specifications.


When we give an input image into a CNN, each of its inner layers generates various activation maps. Activation maps point out the relevant features of the given input image. Each of the CNN neurons generally takes input in the form of a group/patch of the pixel, multiplies their values(colours) by the value of its weights, adds them up, and input them through the respective activation function.

The first (or maybe the bottom) layer of the CNN usually recognizes the various features of the input image such as edges horizontally, vertically, and diagonally.

The output of the first layer is being fed as an input of the next layer, which in turn will extract other complex features of the input image like corners and combinations of edges.

The deeper one moves into the convolutional neural network, the more the layers start detecting various higher-level features such as objects, faces, etc...

Theoretical Experience



We will prepare the project by following the below steps:

- We will be working with Sequential type of modeling, Keras capabilities, image processing techniques
- We will build a web application using the Flask framework.
- Afterwards we will be training our dataset in the IBM cloud and building another model from IBM and we will also test it.

HARDWARE & SOFTWARE Desgining

Hardware Components used

Since we are using the IBM cloud as a platform to execute this project we don't need any hardware components other than our system.

Software Components Used

We will be using Anaconda Navigator which is installed in our system and Watson studio from the IBM cloud to complete the project.

★ **Anaconda Navigator**

Anaconda Navigator is a free and open-source distribution of the Python and R

programming languages for data science and machine learning related applications. It can be installed on Windows, Linux, and macOS. Conda is an open-source, cross-platform, package management system. Anaconda comes with so very nice tools like JupyterLab, Jupyter Notebook, QtConsole, Spyder, Glueviz, Orange, Rstudio, Visual Studio Code. For this project, we will be using Jupiter notebook and spyder.

★ **Watson Studio**

Watson Studio is one of the core services in Cloud Pak for Data as a Service. Watson Studio provides you with the environment and tools to solve your business problems by collaboratively working with data. You can choose the tools you need to analyze and visualize data, to cleanse and shape data, or to build machine learning models. This illustration shows how the architecture of Watson Studio is centered around the project. A project is a workspace where you organize your resources and work with data.

Watson Studio projects fully integrate with the catalogs and deployment spaces:

- Deployment spaces are provided by the Watson Machine Learning service
You can easily move assets between projects and deployment spaces

Experimental Investigations

In this project, we have deployed our training model using CNN on IBM Watson studio and in our local machine. We are deploying 4 types of CNN layers in a sequential manner , starting from

- ★ Convolutional layer 2D:A 2-D convolutional layer applies sliding convolutional filters to 2-D input. The layer convolves the input by moving the filters along the input vertically and horizontally and computing the dot product of the weights and the input, and then adding a bias term.
- ★ Pooling Layer :Pooling layers are used to reduce the dimensions of the feature maps. Thus, it reduces the number of parameters to learn and the amount of computation performed in the network. The pooling layer summarises the features present in a region of the feature map generated by a convolution layer.
- ★ Fully-Connected layer :After extracting features from multiple convolution layers and pooling layers, the fully-connected layer is used to expand the connection of all features. Finally, the SoftMax layer makes a logistic regression classification. Fully-connected layer transfers the weighted sum of the output of the previous layer to the activation function.
- ★ Dropout Layer :There is usually a dropout layer before the fully- connected layer.

The dropout layer will temporarily disconnect some neurons from the network according to the certain probability during the training of the convolution neural network, which reduces the joint adaptability between neuron nodes, reduces overfitting, and enhances the generalization ability.

Project Flow

- User interacts with User interface to upload image
- Uploaded image is analyzed by the model which is integrated
- Once model analyses the uploaded image, the prediction is showcased on the UI

To accomplish this, we have to complete all the activities and tasks listed below.

- Data Collection.
 - Collect the dataset or Create the dataset
- Data Preprocessing.
 - Import the ImageDataGenerator library
 - Configure ImageDataGenerator class
 - Apply ImageDataGenerator functionality to Train dataset and Test dataset
- Model Building
 - Import the model building Libraries
 - Initializing the model
 - Adding Input Layer
 - Adding Hidden Layer
 - Adding Output Layer
 - Configure the Learning Process
 - Training and testing the model
 - Optimize the Model
 - Save the Model
- Application Building
 - Create an HTML file
 - Build Python Code
- Training model on IBM Cloud.

Project Structure

Project folder




data




train and test folders will contain the images to train and test the model respectively.

JupyterNotebook Traning:





IBM Training



Flask

contains the

- html files in templates folder
- .css, .js and image files that are required for html templates in static
- CNN models and
- Python file to run flask



- We are building a Flask Application that needs HTML pages stored in the templates folder and a python script app.py for serverside scripting
- we need the model which is saved and the saved model in this content is ECG.h5
- The static folder will contain js and CSS files.
- Whenever we upload an image to predict, that images are saved in the uploads folder.

WORKING ON PROJECT

1. Dataset Collection

The dataset required for our project "**ECG- Image Based Heartbeat Classification For Arrhythmia Detection Using IBM Watson Studio**" is downloaded from internet.

data:

Name	Date Modified
data	03-07-2022 09:43
test	03-07-2022 09:43
> Left Bundle Branch Block	03-07-2022 09:37
> Normal	03-07-2022 09:39
> Premature Atrial Contraction	03-07-2022 09:40
> Premature Ventricular Contractions	03-07-2022 09:41
> Right Bundle Branch Block	03-07-2022 09:43
> Ventricular Fibrillation	03-07-2022 09:43
train	03-07-2022 10:59
> Left Bundle Branch Block	03-07-2022 10:58
> Normal	03-07-2022 10:58
> Premature Atrial Contraction	03-07-2022 10:59
> Premature Ventricular Contractions	03-07-2022 10:59
> Right Bundle Branch Block	03-07-2022 10:59
> Ventricular Fibrillation	03-07-2022 10:59
.DS_Store	03-07-2022 09:43

2. Image Processing

The dataset of our project includes images of different classes (6 classes) which need "Image Processing" step before feeding the input to ANN. This Image processing include

- Augmenting the image feature, image data generator library.
- Load dataset.
- Apply augmented feature to train set and test set.

3. Model Building

As the Image Processing step is done, now we start to build the CNN model for prediction.

This step includes the following steps:

- Import libs
- Initialize the model
- Add CNN layers
- Configure your learning
- Fit the data
- Save the model

4. Application Building

In this project, we not only build a CNN model, but we also build an Application, a proper platform to use this trained model, using a micro web framework called **FLASK**.

For this we create HTML pages (with associated CSS and JS pages) which shares some information about Arrhythmia and classifications, and in the 'predict' page, we predict the type of Arrhythmia given as Image Input.

These HTML pages are used to make a proper interface (a website) to showcase our project output.

Flask folder includes

templates (to store the .html files needed for our website)

static (to store .css .js and some images included in .html files)

uploads (to store the images we uploaded while using the app)

app.py (flask application development python file)

<modelname>.h5 file (CNN trained model that is given to flask application as input)

5. Train the model on IBM

Name	Type	Size	Authors
📁 .ipynb_checkpoints	File folder		
📄 ECG_IBM.h5	H5 File	10,185 KB	
💻 IBM Model Building.ipynb	IPYNB File	27 KB	
💻 IBM.py	PY File	1 KB	
💻 IBM_Training.ipynb	IPYNB File	14 KB	
RAR my_model.tar.gz	WinRAR archive	6,740 KB	

1. Create a Project using Watson Studio



New project

Define details

Name: CNN_ECG_Sowjanya

Description: ECG- Image Based Heartbeat Classification For Arrhythmia Detection Using IBM Watson Studio

Storage: Cloud Object Storage-tj

Choose project options

Restrict who can be a collaborator ⓘ

Project includes integration with [Cloud Object Storage](#) for storing project assets.

[Cancel](#) [Create](#)

IBM Watson Studio Search in your workspaces Buy Lakshmi Sowjanya Uppana... LU

Projects

Find a project

<input type="checkbox"/>	Name	Date created	Your role	Collaborators	⋮
<input type="checkbox"/>	CNN_ECG_Sowjanya	6 hours ago	Admin	1	⋮

Items per page: 20 1–1 of 1 items 1 1 of 1 pages ⋮

Required IBM services and software

^ Services and software (2)

 Machine Learning-b8	Default	Dallas	Machine Learning	Active	—
 Watson Studio-fy	Default	Dallas	Watson Studio	Active	—

2. Create a model through coding.

IBM Watson Studio

Search in your workspaces

Buy Lakshmi Sowjanya Uppana...

Projects / CNN_ECG_Sowjanya / CNN_ECG_Sowjanya

ECG-Arrhythmia Image Classification

Image Processing

```
In [14]: from tensorflow.keras.models import Sequential
from tensorflow.keras.layers import Dense
from tensorflow.keras.layers import Convolution2D
from tensorflow.keras.layers import MaxPooling2D
from tensorflow.keras.layers import Flatten

In [15]: from tensorflow.keras.preprocessing.image import ImageDataGenerator

In [16]: train_datagen = ImageDataGenerator(rescale = 1./255,shear_range = 0.2,zoom_range = 0.2,horizontal_flip = True)
test_datagen = ImageDataGenerator(rescale = 1./255)

In [19]: #performing data augmentation to train the dataset
x_train=train_datagen.flow_from_directory(directory=r'/home/wsuser/work/data/train',target_size=(64,64),batch_size=32,class_mode='categorical')
#performing augmentation to test the dataset
x_test=test_datagen.flow_from_directory(directory=r'/home/wsuser/work/data/test',target_size=(64,64),batch_size=32,class_mode='categorical')

Found 1534 images belonging to 6 classes.
Found 6825 images belonging to 6 classes.

In [20]: x_train.class_indices
Out[20]: {'left bundle branch block': 0,
'Normal': 1,
'Premature Atrial Contraction': 2,
'Premature Ventricular Contractions': 3,
```

IBM Watson Studio

Search in your workspaces

Buy Lakshmi Sowjanya Uppana...

Projects / CNN_ECG_Sowjanya

Overview Assets Jobs Manage

Find assets

2 assets

All assets

Name	Last modified
CNN_ECG_Sowjanya	3 hours ago Lakshmi Sowjanya Uppanapalli (You)
data.zip	6 hours ago Lakshmi Sowjanya Uppanapalli (You)

Asset types

- > Data 1
- > Source Code 1

Add asset New asset

Data in this project

Drop data files here or browse for files to upload

Items per page: 20 1–2 of 2 items 1 of 1 pages

3. Deploy the model.

The screenshot shows the IBM Watson Studio interface. At the top, there's a navigation bar with 'IBM Watson Studio' on the left, a search bar 'Search in your workspaces' in the center, and user information 'Buy' and 'Lakshmi Sowjanya Uppana...' on the right. Below the navigation bar is a section titled 'Deployments' with a sub-section '1 space'. There are two tabs: 'Activity' and 'Spaces', with 'Spaces' being the active tab. A search bar 'Which deployment space are you looking for?' is present. A table lists the single deployment space: 'ECG_Classification' (Last modified: Jul 5, 2022 6:02 PM, Your role: Admin, Collaborators: 1, Tags: 0, Online deployments: 1, Jobs: 0).

4. Using API key, download the .tar file using Local Jupyter Notebook.

A screenshot of a local Jupyter Notebook file browser. The table lists the following files:

Name	Type	Size	Authors
.ipynb_checkpoints	File folder		
ECG_IBM.h5	H5 File	10,185 KB	
IBM Model Building.ipynb	IPYNB File	27 KB	
IBM.py	PY File	1 KB	
IBM_Training.ipynb	IPYNB File	14 KB	
my_model.tar.gz	WinRAR archive	6,740 KB	

5. Extract .h5 model build in IBM and store it in Flask folder to use it in application created.

A screenshot of a local Jupyter Notebook file browser. The table lists the following files:

Name	Type	Size
static	File folder	
templates	File folder	
uploads	File folder	
.DS_Store	DS_STORE File	9 KB
app.py	PY File	3 KB
ECG.h5	H5 File	9,985 KB
ECG_IBM.h5	H5 File	10,185 KB
ibm_app.py	PY File	3 KB

6. Run the flask application with JupyterNotebook model and IBM model.

```
Anaconda Prompt (anaconda3) - python ibm_app.py


(base) C:\Users\laksh>cd "OneDrive\Desktop\smartinternz courefiles\GUIDED PROJECT\Flask"
(base) C:\Users\laksh\OneDrive\Desktop\smartinternz courefiles\GUIDED PROJECT\Flask>python ibm_app.py
2022-07-05 21:10:20.683270: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cudart64_110.dll'; dlerror: cudart64_110.dll not found
2022-07-05 21:10:20.683567: I tensorflow/stream_executor/cuda/cudart_stub.cc:29] Ignore above cudart dlerror if you do not have a GPU set up on your machine.
2022-07-05 21:10:24.802463: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cudart64_110.dll'; dlerror: cudart64_110.dll not found
2022-07-05 21:10:24.803449: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cublas64_11.dll'; dlerror: cublas64_11.dll not found
2022-07-05 21:10:24.804300: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cublasLt64_11.dll'; dlerror: cublasLt64_11.dll not found
2022-07-05 21:10:24.805149: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cufft64_10.dll'; dlerror: cufft64_10.dll not found
2022-07-05 21:10:24.806046: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'curand64_10.dll'; dlerror: curand64_10.dll not found
2022-07-05 21:10:24.806913: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cusolver64_11.dll'; dlerror: cusolver64_11.dll not found
2022-07-05 21:10:24.807765: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cusparse64_11.dll'; dlerror: cusparse64_11.dll not found
2022-07-05 21:10:24.808609: W tensorflow/stream_executor/platform/default/dso_loader.cc:64] Could not load dynamic library 'cudnn64_8.dll'; dlerror: cudnn64_8.dll not found
2022-07-05 21:10:24.808736: W tensorflow/core/common_runtime/gpu/gpu_device.cc:1850] Cannot dlopen some GPU libraries. Please make sure the missing libraries mentioned above are installed properly if you would like to use GPU. Follow the guide at https://www.tensorflow.org/install/gpu for how to download and setup the required libraries for your platform.



Skipping registering GPU devices...
2022-07-05 21:10:24.809147: I tensorflow/core/platform/cpu_feature_guard.cc:193] This TensorFlow binary is optimized with oneAPI Deep Neural Network Library (oneDNN) to use the following CPU instructions in performance-critical operations: AVX AVX2
To enable them in other operations, rebuild TensorFlow with the appropriate compiler flags.
* Serving Flask app "ibm_app" (lazy loading)
```

OUTPUT of Flask App Development

about.html

According to the World Health Organization (WHO), cardiovascular diseases (CVDs) are the number one cause of death today. Over 17.7 million people died from CVDs in the year 2017 all over the world which is about 31% of all deaths, and over 75% of these deaths occur in low and middle income countries. Arrhythmia is a representative type of CVD that refers to any irregular change from the normal heart rhythms. There are several types of arrhythmia including

- *Ventricular Fibrillation*

- *Premature Atrial Contraction*

- *Premature Ventricular Contraction*

- *Right Bundle Branch Block*

- *Left Bundle Branch Block*


• **Left Bundle Branch Block**

Although single arrhythmia heartbeat may not have a serious impact on life, continuous arrhythmia beats can result in fatal circumstances. Electrocardiogram (ECG) is a non-invasive medical tool that displays the rhythm and status of the heart. Therefore, automatic detection of irregular heart rhythms from ECG signals is a significant task in the field of cardiology.

The prevalence and mortality rates of cardiovascular disease (CVD) continue to rise. As a result, frequent cardiac rhythm monitoring has become an increasingly critical and vital aspect of managing and preventing CVDs. The automatic diagnosis of cardiac illness relies heavily on the classification of electrocardiogram signals. A stroke can result in brain damage and necessitates immediate medical attention. To diagnose an arrhythmia, a doctor must first recognize the abnormal heartbeat and attempt to determine its cause or trigger. Thanks to the development of artificial intelligence and Science that has enabled us to predict the cases of arrhythmia far better than doctors by the use of Convolutional Neural Networks. We in this project aim to diagnose the type of arrhythmias by the spectrograms of numerical data of the ECG images.

info.html

ECG

Electrocardiogram

NORMAL

Note that heart is beating in a regular sinus rhythm between 60-100 beats per minute (specifically 82bpm).

Normal Sinus Rhythm

Information about Arrhythmia x +

localhost:5000/info

Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...


NORMAL

Note that heart is beating in a regular sinus rhythm between 60-100 beats per minute (specifically 82bpm).

All important intervals on this recording are within normal ranges.

The normal ECG patterns seen in children differ considerably from those in adults.

Normal Sinus Rhythm



Adobe Stock #33978108

Information about Arrhythmia x +

localhost:5000/info

Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...


VENTRICULAR FIBRILLATION

- A life-threatening heart rhythm that results in a rapid, inadequate heartbeat.
- Ventricular fibrillation (VF) is a rapid, life-threatening heart rhythm starting in the bottom chambers of the heart. It can be triggered by a heart attack.
- Because the heart doesn't pump adequately during ventricular fibrillation, sustained VF can cause low blood pressure, loss of consciousness or death.
- Emergency treatment includes immediate defibrillation with an automated external defibrillator (AED) and cardiopulmonary resuscitation (CPR). Long-term therapy includes implantable defibrillators and medications to prevent recurrence.

Normal Sinus Rhythm Steady rhythm and rate



Ventricular Fibrillation Chaotic rhythm and rate



©2021 Cleveland Clinic

Information about Arrhythmia x +

localhost:5000/info


Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...

PREMATURE ATRIAL CONTRACTION

Premature atrial contractions (PACs) are extra heartbeats that start in the upper chambers of your heart. When the premature, or early, signal tells the heart to contract, there may not be much blood in the heart at that moment. That means there's not much blood to pump out.

SYMPTOMS:

- You may not have any symptoms. If you do, your symptoms may include heartbeats that:
- Have a lot of force now and then.
- Skip.
- Pound (palpitations).
- You may also have anxiety or shortness of breath.



An ECG strip showing multiple leads (I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6). The strip displays a regular rhythm with occasional提前的搏动 (premature beats) that appear earlier than expected. These premature beats are characterized by a narrow QRS complex followed by a normal T wave.

Information about Arrhythmia x +

localhost:5000/info


Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...

PREMATURE VENTRICULAR CONTRACTION

Premature ventricular contractions (PVCs) are extra heartbeats that begin in one of the heart's two lower pumping chambers (ventricles). These extra beats disrupt the regular heart rhythm, sometimes causing a sensation of a fluttering or a skipped beat in the chest.

SYMPTOMS:

- Fluttering.
- Pounding or jumping.
- Skipped beats or missed beats.
- Increased awareness of the heartbeat.



An ECG strip showing multiple leads (I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6). The strip shows a regular rhythm with a single premature ventricular contraction (PVC). This is identified by a wide, distorted QRS complex followed by a normal T wave. The subsequent beats return to the regular sinus rhythm.

Information about Arrhythmia ▾ +

localhost:5000/info

Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...

RIGHT BUNDLE BRANCH BLOCK

Right bundle branch block is an obstacle in your right bundle branch that makes your heartbeat signal late and out of sync with the left bundle branch, creating an irregular heartbeat.

SYMPTOMS:

In most people, bundle branch block doesn't cause symptoms. Some people with the condition don't know they have bundle branch block.

Rarely, symptoms of bundle branch block may include fainting (syncope) or feeling as if you're going to faint (presyncope).

When to see a doctor?

If you've fainted, see a health care provider to rule out serious causes.

If you have heart disease or have been diagnosed with bundle branch block, ask your provider how often you should have follow-up visits.



Information about Arrhythmia ▾ +

localhost:5000/info

Microsoft account [...] Studies VIT APPS My Drive vit mail vtop Home Page - Select...

LEFT BUNDLE BRANCH BLOCK


A left bundle branch block usually is a sign of an underlying heart disease, including dilated cardiomyopathy, hypertrophic cardiomyopathy, high blood pressure, aortic valve disease, coronary artery disease and other heart conditions. While left bundle branch block can appear in healthy people, it most often does not.

SYMPTOMS:


- Left bundle branch block often doesn't cause any symptoms. fatigue and shortness of breath can be considered as symptoms for it.
- A delay in the arrival of electrical impulses to the heart's left ventricle can cause syncope (fainting), due to unusual heart rhythms that affect blood pressure.
- Some people might also experience something called presyncope. This involves feeling like you're about to faint, but never actually fainting.




upload.html



CNN Prediction:





Advantages, Disadvantages and Applications

Advantages

- The proposed model predicts Arrhythmia in images with a high accuracy rate of nearly 96%
- The early detection of Arrhythmia gives better understanding of disease causes, initiates therapeutic interventions and enables developing appropriate treatments.

Disadvantages

- Not useful for identifying the different stages of Arrhythmia disease.
- Not useful in monitoring motor symptoms

Applications

- It is useful for identifying the arrhythmia disease at an early stage.
- It is useful in detecting cardiovascular disorders

Conclusion

- Cardiovascular disease is a major health problem in today's world.
- The early diagnosis of cardiac arrhythmia highly relies on the ECG.
- Unfortunately, the expert level of medical resources is rare, visually identify the ECG signal is challenging and time-consuming.
 - The advantages of the proposed CNN network have been put to evidence.
 - It is endowed with an ability to effectively process the non-filtered dataset with its potential anti-noise features. Besides that, ten-fold cross-validation is implemented in this work to further demonstrate the robustness of the network.

Future Scope

For future work, it would be interesting to explore the use of optimization techniques to find a feasible design and solution. The limitation of our study is that we have yet to apply any optimization techniques to optimize the model parameters and we believe that with the implementation of the optimization, it will be able to further elevate the performance of the proposed solution to the next level.

References

- <https://www.analyticsvidhya.com/blog/2021/05/convolutional-neural-networks-cnn/>
- <https://www.mathworks.com/help/deeplearning/ref/nnet.cnn.layer.convolution2dlayer.html?jsessionid=0a7e3bc26fabda07a5032030294b>
- https://smartinternz.com/externship_dyn/1/artificial-intelligence

Project Links

👉 **Github Link :**

<https://github.com/smartinternz02/SI-GuidedProject-78150-1656680500/tree/main>

👉 **Demo Link:**

https://drive.google.com/file/d/1VaZoEufjwGuQ2PuOTCzTQepT_cY6i52R/view?usp=sharing

THE END
